

**EASY
PRECISE
VERSATILE
FAST**

primosplint




light cured composite splint material

in general light cured composite is the most precise,
distortion-free and antiallergenic material for
night guards and more

EASY PRECISE VERSATILE FAST



light cured
composite
rods, MMA
and peroxide
free





primosplint

[uses]

- function therapy bite splints
- night guards
- verification indexes
- insertion guides
- transfer jigs

[advantages]

- very quick and easy to work with
- modelling playdough like consistency at room temperature – no heat needed
- cures reliably in most conventional light curing units
- extremely precise, no noticeable shrinkage, expansion or distortion
- abrasion-optimized, protects the tooth structure and the TMJ
- MMA and peroxide free – lowest allergic potential for technician and patient
- fully autoclavable – no distortion during steam sterilisation
- can be easily bonded to other materials (PMMA-acrylic or suck down matrix)
- no taste or smell and easy to repair

technical tips

When opening a new stick of Primosplint gently knead it for about 10 to 15 seconds to soften up the material. This will make the material softer and easier to work with.

For maximum clarity lightly sandblast the finished splint (approx. 110 micron aluminium oxide, 30 PSI) inside and outside. Then apply one thin layer of primoglaze on the inside, blow off any excess and apply two thin layers on the outside. Light cure according to the instructions.

If adding Primosplint to a suck down matrix, make sure that the surface is roughened well enough before applying Primostick bonder.

Do not use strong strobe lights for curing.

Any further technical questions – please call **203-682-6429**

background

Primosplint is a thixotropic material, meaning the more energy put into it, the softer and more pliable it gets.

Primosplint provides the highest strength without any brittleness because tensile strength and modulus of elasticity have been optimized for bite splint applications.

Splints made of primosplint have a very high accuracy of fit, meaning they are utmost dimensionally stable during and after polymerization. This prevents "rocking" or "seesawing" effects.

The mechanical properties of primosplint have been carefully tuned to relieve the temporomandibular joint and protect the tooth structure. With a flexural strength of > 70 MPa and a modulus of elasticity of > 1500 MPa, bite splints made of primosplint meet the mechanical requirements asked for by modern splint therapy.



www.primotecusa.com

JUST A FEW STEPS TO PERFECT RESULTS!!!



Finished Case

[working procedure]



1
primosplint - light cured composite in rod shape



2
easy adaptation at room temperature (no heat needed)



3
establish centric and disclusion



4
polymerisation in one of the metalight units...



...or in most conventional UV-A curing lights



5
lift-off the model



6
no shrinkage, no rocking, easy to finish and glaze



7
verify occlusion and glaze or polish

ORDERING INFO

primosplint Starter Kit

Item #	Item	Price
PS-KIT	Primosplint PS-Kit1 – one each of: PS 100, box with 12 primosplint rods PS 202, bottle of light cure bonding material PS 203, bottle of stone model separator PS 204, bottle of light cure glaze lacquer PS 205, bottle of cleaning liquid IM 200, instruments and finger separator	\$ 249.00



primosplint refill components


Item #	Item	Price
PS 100	Box with 12 rods of Primosplint – light cure composite material – 8.3 grams each per pack	\$ 129.00
PS 202	10 ml bottle primostick – light cure bonding material	\$ 30.00
PS 203	50 ml bottle primosep – stone model separator	\$ 27.00
PS 204	10 ml bottle primoglaze – light cure glaze lacquer	\$ 30.00
PS 205	50 ml bottle primoclean – cleaning liquid	\$ 27.00
IM 200	Metatouch – Instruments and fingers separator, 43g jar	\$ 45.00



publications

- **Implant Cases Made Easy**
by Christina Ketzinger CDT
and Andreas Hoffmann MDT
published in Spectrum dialogue
Vol. 7 No. 8
- **primosplint: for faster, easier splints**
by Andreas Hoffmann MDT
published in Dental Technologies No. 71

Find these and more publications as well as the Material Safety Data Sheets (MSDS) under www.primotecusa.com or call primotec at **203-582-6429** for your free printed copy.



Dr. Diether Reusch

Privatpraxis für Zahnmedizin,
Westerburg, Germany

"Especially for our CMD therapy patients, primosplint offers us the needed perfect passive fit to ensure best treatment results. To easily verify occlusal contacts and disclusion we prefer the primosplint surface finish with lesser translucency."



[easy, precise, versatile, fast]

light cured composite rods, MMA and peroxide free –
for function therapy bite splints and many more applications.



Primotec USA
215 Westport Ave.
Norwalk, CT 06851
203-682-6429

info@primotecusa.com
www.primotecusa.com