

die stone

DRS 100

- stays firm on the spoon and flows nicely on the vibrator due to its high thixotropy and stability
- shatterproof sawing, trimming and drilling for ultra-smooth and sharp model edges

resin-reinforced, super hard type 4 die stone

primostone DRS 100



- **Box 55 Lb.**
- **US\$ 129.00**

primostone DRS 100 is the uniquely formulated, high end die stone from primotec. It is resin-reinforced and highly thixotropic for fine, smooth dies as well as opposing, working and study models.

- **available in gold-brown or vanilla**

die stone

DCL 50

- very thixotropic - stays on the spoon and flows on the vibrator
- removal from the mold after 30 minutes due to fast setting time

classic high-performance type 4 die stone

primostone DCL 50



- **Box 55 Lb.**
- **US\$ 109.00**

primostone DCL 50 is the classic super hard die stone from primotec. It is highly abrasion resistant and super precise due to its very high final hardness. Suitable for all dies, opposing, working and study models.

- **available in gold-brown or ivory**

implant stone

IXO

- no dimensional change of the model due to 0% setting expansion
- 100% precision for implant and long-span cases

ZERO-expansion, super hard type 4 implant model stone

primostone IXO



- **Box 50 Lb.**
- **US\$ 119.00**

primostone IXO is the world's first proven 0,00% expansion super hard type 4 stone from primotec. It is ideal to use for implant case models and all applications that require maximum long-span precision. primostone IXO is NOT used for dies.

- **available in gold-brown or ivory**

base stone

MBC 100

- precise interproximal contacts and best fit due to very low expansion of the model base
- extra smooth denture surfaces when used for flasking procedures

low expansion, super hard type 4 base stone

primostone MBC 100

primostone MBC 100 is the special, readily flowing type 4 base stone from primotec. It is mainly used for basing crown & bridge models, but also works well for split cast models, full or partial denture models and denture flasking procedures.

- available in royal-blue or white



lab stone

FCM 100

- smooth, creamy mixture for accurate bubble-free models
- heat resistant, not effected by boil-out procedures

versatile, all-purpose hard type 3 lab stone

primostone FMC 100

primostone FCM 100 is the hard, general-purpose type 3 lab stone from primotec. It is best to use for study and working models, for full and partial denture models, denture flasking procedures and partial repairs.

- available in blue



arti plaster

AOX

- for the most precise transfer of occlusal relation from patient to articulator
- excellent adhesion to the mounted models

ZERO-expansion articulator mounting plaster

primostone AOX

the world's first and only 0.00% expansion articulator mounting plaster from primotec. This makes AOX 100 the perfect articulating plaster, since it does not cause any bite increase and maintains optimal occlusion regardless of how much plaster is used.

- available in white

