

# PREMIO Milling burs for Sirona®

## PREMIO ZR DIAMOND

Zirconia



High-End Milling burs for zirconia

- with **long life diamond coating** for maximum longevity
- 3 flutes
- 3 mm shank Ø

DIAMOND

Product #	Name	Price
<b>SI.3.05.D</b>	0,5 mm Premio ZR Diamond (2 flutes)	<b>\$119.00</b>
<b>SI.3.10.D</b>	1,0 mm Premio ZR Diamond	<b>\$119.00</b>
<b>SI.3.25.D</b>	2,5 mm Premio ZR Diamond	<b>\$119.00</b>

Zirconia  
PMMA  
Wax



Standard burs for zirconia, PMMA and wax

- **uncoated** long life
- 2 or 4 flutes
- 3 mm shank Ø

STANDARD

Product #	Name	Price
<b>SI.3.05</b>	0,5 mm Premio ZR Standard	<b>\$49.00</b>
<b>SI.3.10.ZR</b>	1,0 mm Premio ZR Standard	<b>\$49.00</b>
<b>SI.3.25.ZR</b>	2,5 mm Premio ZR Standard (4 flutes)	<b>\$49.00</b>

PMMA  
Wax



Standard burs for PMMA, wax and composite

- **uncoated** long life
- 1 or 2 flutes
- 3 mm shank Ø

STANDARD

Product #	Name	Price
<b>SI.3.05</b>	0,5 mm Premio PM Standard (2 flutes)	<b>\$49.00</b>
<b>SI.3.10</b>	1,0 mm Premio PM Standard	<b>\$49.00</b>
<b>SI.3.25</b>	2,5 mm Premio PM Standard	<b>\$49.00</b>

# PREMIO Milling burs for Sirona®

## PREMIO HC STANDARD

Ceramic  
Blocks®



Precise milling burs for ceramic and hybrid-ceramic blocks

- with **diamond coating** for outstanding longevity
- galvanically bonded
- shank with special carbide metal core 3 mm Ø

# STANDARD

Product #	Name	Price
<b>SI.3.12.EMA</b>	1,2 mm Premio HC Standard	<b>\$39.00</b>
<b>SI.3.14.EMA</b>	1,4 mm Premio HC Standard	<b>\$39.00</b>
<b>SI.3.22.EMA</b>	2,2 mm Premio HC Standard	<b>\$39.00</b>