





PRICE LIST  
**PREMIOtemp PMMA Blanks**



		<b>MULTICOLOR</b>	
		<p><b>High-End PMMA Blanks</b></p> <p><b>MULTI</b></p> <p>For highly aesthetic long term temporary restorations</p>	<ul style="list-style-type: none"> <li>• Impressive aesthetics by 5 layers of color</li> <li>• Multilayered with no separation between layers</li> <li>• Natural fluorescence</li> <li>• Also available in AG-Shape</li> </ul>
Diameters:	98mm with 2 steps 95mm (ZZ-Shape) 71mm (AG-Shape)		
Heights:	16mm, 20mm or 25mm (30mm on request)		
Colors:	A1 to D4, plus Bleach (BL1, BL2, BL3, BL4) Pink (Gingiva) and Bicolor		
<b>Price:</b>	<b>\$99.00</b> (16mm & 20mm) <b>\$109.00</b> (25mm) <b>\$119.00</b> (30mm)		

		<b>MONOCOLOR</b>	
		<p><b>High-End PMMA Blanks</b></p> <p><b>MONO</b></p> <p>Good aesthetics for regular short term provisionals</p>	<ul style="list-style-type: none"> <li>• Excellent combination of color and translucency</li> <li>• Easy to mill</li> <li>• Highly biocompatible</li> </ul>
Diameters:	98mm with 2 steps 95mm (ZZ-Shape)		
Heights:	16mm, 20mm or 25mm (30mm on request)		
Colors:	A1 to D4, plus Bleach (BL1, BL2, BL3, BL4) and Pink (Gingiva)		
<b>Price:</b>	<b>\$49.00</b> (18mm) <b>\$59.00</b> (25mm) <b>\$69.00</b> (30mm)		

		<b>CLEAR</b>	
		<p><b>High-End PMMA Blanks</b></p> <p><b>CLEAR</b></p> <p>State of the art for milled bite splints</p>	<ul style="list-style-type: none"> <li>• Glossy surface directly after milling</li> <li>• Transparent as glass</li> <li>• High strength (&gt;130MPa)</li> </ul>
Diameters:	98mm with 2 steps 95mm (ZZ-Shape)		
Heights:	10, 12, 14, 16, 18, 20, 25mm		
Colors:	Transparent		
<b>Price:</b>	<b>\$29.00</b> (98mm x 10mm-20mm) <b>\$35.00</b> (98mm x 25mm)		
			<p>Certified as Class 2A medical devices according to EC Directive 93/42/EEC</p>

PRICE LIST  
**PREMIOtemp PMMA Blanks**



**CLEAR FLEX**



**High-End PMMA Blanks**

**CLEAR FLEX**

The thermoplastic material for milled bite splints

- Innovative milling blank
- Automatically adapts to the patient's dental situation when warmed up
- Tension-free wearing comfort
- Durable due to maximum break resistance

Diameters: 98mm with 2 steps  
 95mm (ZZ-Shape)

Heights: 16mm, 20mm or 25mm

Colors: Translucent

**Price:** **\$69.00**

Certified as Class 1

**ORDER**

- by Email: [info@primotecusa.com](mailto:info@primotecusa.com)
- by fax: 212-929-7444
- or by phone: 203-682-6429

**Your Order:**

**1** Type **MULTI**

Type **MONO**

Type **CLEAR**

Type **CLEAR FLEX**

**2** Diameter: \_\_\_\_\_

Height: \_\_\_\_\_

Color: \_\_\_\_\_

Quantity: \_\_\_\_\_

If you are already a **primotec customer**, please insert only your **name** and **customer number**:

Name / Lab: \_\_\_\_\_ Customer number: \_\_\_\_\_

**Please send us our order to following address\*:**

Name / Lab: \_\_\_\_\_

Contact person / Position: \_\_\_\_\_

Street: \_\_\_\_\_

ZIP / City: \_\_\_\_\_

Telephone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Date / Signature: \_\_\_\_\_

**\*Privacy:**

At any time you have a right for information, rectification and cancellation regarding the processing and use of your data. For sure, you can revoke your consent at any time at primotec.