


# PREMIO Milling burs for **Wieland**® 5 axis-systems

## PREMIO ZR DIAMOND

Zirconia



High-End Milling burs for zirconia

- with **long life diamond coating** for maximum longevity
- 3 flutes
- 3 mm shank Ø

DIAMOND

Product #	Name	Price
<b>VHF.3.10.D</b>	1,0 mm Premio ZR Diamond	<b>\$119.00</b>
<b>VHF.3.25.D</b>	2,5 mm Premio ZR Diamond	<b>\$119.00</b>

## PREMIO PM BLUE

PMMA  
PEEK  
Zirconia



High-End Milling burs for PMMA, PEEK and zirconia

- with **Durablue coating** for outstanding longevity
- 2 flutes
- 3 mm shank Ø

BLUE

Product #	Name	Price
<b>VHF.3.10.B5</b>	1,0 mm Premio PM Blue	<b>\$59.00</b>
<b>VHF.3.25.B5</b>	2,5 mm Premio PM Blue	<b>\$59.00</b>

## PREMIO PM STANDARD

PMMA  
PEEK  
Wax



Standard burs for PMMA, PEEK and wax

- **uncoated** long life
- 2 flutes
- 3 mm shank Ø

STANDARD

Product #	Name	Price
<b>VHF.3.07.5</b>	0,7 mm Premio PM Standard	<b>\$49.00</b>
<b>VHF.3.10.5</b>	1,0 mm Premio PM Standard	<b>\$49.00</b>
<b>VHF.3.25.5</b>	2,5 mm Premio PM Standard	<b>\$49.00</b>

# PREMIO Milling burs for **Wieland**<sup>®</sup> 5 axis-systems

## PREMIO PM STANDARD

PMMA  
PEEK  
Wax



Standard burs for PMMA, PEEK and wax

- **uncoated** long life
- 1 flute
- 3 mm shank Ø

# STANDARD

Product #	Name	Price
<b>VHF3.07.PM5</b>	0,7 mm Premio PM Standard	<b>\$49.00</b>
<b>VHF3.10.PM5</b>	1,0 mm Premio PM Standard	<b>\$49.00</b>
<b>VHF3.25.PM5</b>	2,5 mm Premio PM Standard	<b>\$49.00</b>